


Data sheet

Hospital/ Emergencys/ Ambulance/ Nursing home
Disposable Transfer Sheet Multi




Product Code	193629DK	
Manufacturer	EMS MEDICAL A/S	
Made in	Denmark	
Product Description	Absorbent transfer sheet with 8 integrated grip holes for easy and comfort lifting/ transfer the patient up to 130 KG.. Waterproof, and good breathable	
Features	With 8 integrated grip holes	
Size	Outer measure: 80 x 190 +/-3 cm	Inner core Absorbent: 65 x 170 +/-3 cm
Packaging	1 Pcs/ Bag,	
	30 Pcs/ Ctn.	
	Carton size 58cm*35cm*50cm	
	Gross weight 9.2kg, Net weight 10,8kg.	
Shelf life	5 years	
Storage	Keep product in-House, Protect from sunlight, dust and humidity, temperature: 5° to 40° C	
Usage	Single use	
Sterilization	Non Sterile	Ready for EO sterile in Hospitals/Clinics before use if needed. Packaged with color indicator for sterilization.
	Product classification Medical device acc. to MDR 2017/745	
Material	Top layer	Non-woven 20G/m ²
	Middle Layer	Tissue 13G/m ² + Fluff 24G/m ² + Tissue 13G/m ²
	Bottom Layer	SS Non-woven 105G/m ² , with PE film 18G/m ²
	Total	380G/m ² , 8G/m ² SAP
	Colour	Blue / White
Absorbency	2,3 L fluid uptake	
Classify	Medical Device, Class I	
Free of	Glass fibres, latex, additives	
Max load for transfer	130kg.	
Languages	Danish, English, Swedish, Norwegian, French, German, Spanish, Portuguese, etc.	



Liquid Absorbency



Load 130KG.



Extra Durable



For further information:

EMS MEDICAL A/S, DENMARK

+45-52714444

www.emsmedical.dk