

# Datasheet

## Medical Disposable Hair Cover



MDR 2017/745, CLASS 1




<b>Product No.</b>	<b>197122DK</b>	
<b>Manufacturer</b>	EMS MEDICAL A/S	
<b>Made in</b>	Denmark	
<b>Product Description</b>	Medical Disposable Hair Cover	
<b>Features</b>	1x3mm Elastic	
<b>Packaging</b>	1 Euro Pallet = TBC	
<b>Shelf life</b>	5 years	
<b>Storage</b>	Keep product in-House, Protect from sunlight, dust and humidity, temperature: 05° to 40° C	
<b>Usage</b>	Single use	
<b>Sterilization</b>	Usage	Non Sterile
<b>Product classification</b>	Medical device acc. to MDR 2017/745	
<b>Material</b>	Material	Allergy free Spunbund non-woven polypropylene, Extra soft
	Weight/m2	15g
	Colour	White
<b>MD tensile strength</b> Test method NWSP 110.4	Avg	44,8 N
	Min	42,9 N
	Max	46,6 N
<b>CD tensile strength</b>	Avg	20,2 N
	Min	19,9 N
	Max	21,2 N
<b>Free of</b>	Glass fibres, latex, additives	
<b>Performance</b>	Air Permeability (Pa/cm2)	< 5
	Microbial Cleanliness (cfu/g)	< 10
<b>Languages</b>	Danish, English, Swedish, Norwegian, French, German, Spanish, Portuguese, Finnish Other EU languages can be added as a label.	
<b>Material Standard</b>	Asthma-Allergy Friendly	
	Oeko-Tex 100	
	Extra Soft for skin	
	Environmentally Friendly	



Asthma-Allergy Friendly



Sustainability



Extra Durable



Extra Soft



Environmentally Friendly



For further information:

EMS MEDICAL A/S, Jegstrupvej 40, 8361 Hasselager, Denmark

+45-52714444

www.emsmedical.dk