



**Ref No.** 201164DK

**Disposable protective sleeve, with cuffs, Non-woven (SMS)**

- ✓ **With elastic in the upper arm**
- ✓ **Minimize cross contamination**
- ✓ **Cuffs with ribbing**

**STERILE EO**



## Product description

Disposable protective sleeve with cuffs and elastic in the upper arm.  
For use in medical and hygienic environments.

## Specifications

<b>Product name</b>	Disposable protective sleeve
<b>Trademark</b>	EMS Medical
<b>Color</b>	Blue
<b>Material</b>	Non Woven (SMS)
<b>Disposable/ Reusable</b>	Disposable
<b>Thickness</b>	40 gr.
<b>Sizes</b>	One size
<b>Length</b>	64 cm
<b>Features</b>	Sterile, Elastic in the upper arm, Cuffs with ribbing
<b>Certificates</b>	CE
<b>Directives, legislation and regulations</b>	EU 2017/745, 93/ 42/ EEC
<b>Sterilization method</b>	EO Sterile
<b>Classification</b>	Class Is
<b>Storage</b>	Store dry, at room temperature and out of direct sunlight.
<b>Durability</b>	5 years from date of production
<b>Safety regulations and warnings</b>	Material selected with care for minimal risk of allergy and hypersensitivity, but cannot be guaranteed for everyone.
<b>Disposal of product</b>	Can be disposed of with normal household waste. In case contamination is suspected, discard as clinical waste.
<b>Disposal of packaging</b>	Can be recycled or incinerated.
<b>Packaging</b>	PE/PAPER, Flat Sterile pack, Easy to open, Compact size for easy storage
<b>Environment</b>	- Carton 70% recycle / Transport carton 70% Recycle - FSC approved packaging - Sustainability

## Packaging

Unit	Contains	REF No			
Carton	400 pcs	201165			
Sterile Pack	1 pcs				
Pcs	1 pcs				



**STERILE EO**

