



## Disposable Hair cap Non Woven, 15 gr.

- ✓ For use in medical and hygienic environments
- ✓ Medical devices, Class I
- ✓ Soft and comfortable



### Product description

Disposable hair cap's are used to protect the environment, where a high hygienic standard is required.

### Specifications

<b>Product name</b>	Disposable hair caps
<b>Trademark</b>	EMS Medical
<b>Color</b>	White
<b>Material</b>	Polypropylen Non-woven
<b>Disposable/ Reusable</b>	Disposable
<b>Thickness</b>	15 gr.
<b>Size</b>	One size
<b>Classification</b>	Class I
<b>Features</b>	<ul style="list-style-type: none"> <li>- With elastic</li> <li>- Asthma - Allergy friendly</li> <li>- Extra soft</li> </ul>
<b>Certificates</b>	CE for Personal protective equipments (PPE)
<b>Directives, legislation and regulations</b>	MDR 2017/745, (EF) Nr. 10/2011, (EU) 2016/425 (EF) nr. 1935/2004, (EF) Nr.-2023/2006, REG nr. 681 af 25/05/2020
<b>Storage</b>	Store dry, at room temperature and out of direct sunlight.
<b>Durability</b>	5 years from date of production
<b>Safety regulations and warnings</b>	Non Sterile Material selected with care for minimal risk of allergy and hypersensitivity, but cannot be guaranteed for everyone.
<b>Disposal of product</b>	Can be disposed of with ordinary household waste, sorted according to local regulations..
<b>Disposal of packaging</b>	Can be recycled or incinerated
<b>Environment</b>	<ul style="list-style-type: none"> <li>- Carton 70% recycle / Transport carton 70% Recycle</li> <li>- LDPE minimum 80% recycling</li> <li>- FSC approved packaging</li> <li>- Sustainable</li> </ul>



**Packaging**

Pcs/ Bag	Bag/ Carton (Ctn)	Ctn/ Palle	Size	EAN No.	REF No.
100 pcs	15 bag	20 Ctn	One size	5700002197122	197122DK
2000 pcs	1 bag	20 Ctn	One size	5700001801129	801129DK